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Vitrification solution used for mouse embryo cryopreservatio (Mouse-specific)

Product	Item No.	specification	Application
EFS40	72023	1mL	Embryo cryoprotectant solution in mice
EFS20	72024	1mL	Sperm cryoprotectant solution in mice

Composition

DPBS、BSA、Ficoll.

Uses and Principles

The combination of EFS40 and EFS20 can be used for cryopreservation of mouse embryos.

Storage conditions and expiration date

1. Storage Temperature: -20°C, Protect from light.
2. Expiration date: 6 months, unused reagents after opening should be sealed with parafilm and Stored at -20°C.

Embryo selection:

Select embryos with normal morphology for cryopreservation.

Methods

1. Take EFS20 and EFS40 out from -20°C been placed in a metal temperature controller (0°C).
2. Take out the cryovials and place it in a metal temperature controller (0°C).
3. Pipette 50 μ L of EFS40 into the bottom of the cryovials gently.
4. Make one drop in a 35mm dish with EFS40 and several 50 μ L drops with EFS20.
5. 2- cell stage embryos were collected that obtained from IVF.
6. The embryos were placed in drop of EFS20 for 2 min. It was observed that the embryos were suspended in EFS20.
7. At 80s, the pipette can be loaded with a few of EFS40, and followed by aspirating all embryos placed in the EFS20 and gently blowing into the bottom of the cryotube .you can be determined by the entry of bubbles to judge whether all embryos have entered

the cryovial, put it in a metal bath (0 °C) for 1min. Immediately immerse the cryovials with embryos into LN₂ using a forceps.

Cautions:

1. When embryo cryopreservation, do not cover the cryovials cap too tightly to avoid wasting time when the cap cannot be unscrewed during embryo thawing.
2. Ice box can be used instead of metal bath to ensure 0°Ctemperature.
3. Aspirate embryos with a minimal volume of EFS20 medium.
4. Make sure that the cryotube is immersed in the LN₂.
5. It is recommended that embryos be frozen at 50-100R/ cryovials.