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Pregnant Mare Serum Gonadotropin (PMSG) and Human chorionic gonadotropin(HCG) (Mouse-specific)

Product	Item No.	specification	Application
PMSG	72122	1mL	Induce follicular development
HCG	72120	1mL	Induce the rupture of the matured follicles

Composition

PMSG Freeze-dried powder, HCG Freeze-dried powder, Saline

Storage conditions and expiration date

1. Storage: Stored away from light, at -20 °C.
2. Expiration date: 3 months

Mouse

Female mice 3.5-4 weeks of age, [Note]: Efficient induction of superovulation in mice depends on several variables:the age and the does of gonadotropins,and the mice strains used (Table 1).

Instructions

1. After the hormone is completely thawed, gently mix well and draw the PMSG (or HCG) with a disposable 1ml syringe . Evacuating the air from the syringe.
2. Pick up the mouse by the scruff of its neck as close t the ears as possible.Pierce the skin and abdominal muscles wit the hypodermic needle, and inject the solution into the i.p cavity,taking care to aviod the diaphragm and other organs.
3. For most strains,a 46-to 48h interval between PMSG and the hCG injections has been found to be optimal for oocyte yield. Ovulation takes place between 10-13h after injection of hCG, Collect oocyte-cumulus complexes by tearing the oviduct ampulla under a stereoscopic after 13-16 hours of HCG. beyond 16 hours

could result in oocyte aging, which could lead decrease the oocyte fertilization ability, therefore the oviducts are collected around 14 hours after HCG injection.

Cautions:

1. Repeated freezing and thawing of hormones will affect the superovulation effect.
2. Hormone dose selection for different strains and week's old of mouse .

Strains	Age	PMSG	HCG	Notes
C67BL/6J	3.5-4W	5 IU	5 IU	
C67BL/6N	3.5-4W	5 IU	5 IU	
FVB	3.5-4W	5 IU	5 IU	
ICR	3.5-4W	5 IU	5 IU	
Balb/c	8-10W	5 IU	5 IU	Pre-injection of GnRH
NCG/NOD	8-10W	10 IU	10 IU	
SD Rat	4-5W	50IU	50IU	

3. At intraperitoneal injection, stay 2 seconds to confirm that there is no leakage when the mouse were injected ,If leakage occurs, an appropriate amount of injection should be added and then put back into the cage.

4. If the customer need hormone concentration storage, you can contact us separately.